

## Protocol of DNA barcode duckweeds by *atpF-atpH* marker

1. Order primer pairs for marker *atpF-atpH* and prepare them into 5pmol work solution.

Forward: ACTCGCACACACTCCCTTTCC

Reverse: GCTTTTATGGAAGCTTTAACAAT

2. Extract total genomic DNA using CTAB method (see our protocol).

3. Quantify DNA concentration by Nanodrop.

4. PCR mix:

DNA	1ul (50-100ng)
Primer	(1+1)ul (5pmol each)
H <sub>2</sub> O	9.5ul
Red taq mixture	12.5ul
Total	25ul

5. PCR program

94°C 2min

94°C 15s

50°C 15s 35 cycles

72°C 40s

72°C 5min

6. Purify PCR products with ExoSap-IT™.

7. Directly sequence purified PCR products on an ABI3730 automated sequencer using the same primers as in the PCR reactions.

## Chemical

JumpStart™ Redtop® Ready-Mix™ Reaction Mix (P1107, Sigma)

ExoSap-IT™ (USB Corp.)

## Reference

1. Hollingsworth PM FL, Spouge JL, Hajibabaei M, Ratnasingham S, et al (2009) A DNA barcode for land plants. *Proc Natl Acad Sci U S A* 106: 12794-12797.
2. Wang W, Wu Y, Yan Y, Ermakova M, Kerstetter R, et al. (2010) DNA barcoding of the Lemnaceae, a family of aquatic monocots. *BMC Plant Biology* 10: 205.